

Instructions for use of medical tweezers from Scholz Medical GmbH

Manufacturer

Scholz Medical GmbH
 Robert-Bosch-Strasse 1
 83052 Bruckmühl
 Germany
<https://www.scholz-instruments.com/>
info@scholz-medical.de
 +49 8062 8070801



This instruction manual applies to the following products

Description	Product no.	Basic UDI-DI	Category	Rule
Splinter tweezers, 110mm	TABD9649	426052621TABD01X6	Ila	7
Anatomical tweezers, 300mm	NUSIBD2299	426052621TABD01X6	Ila	7
Dental tweezers, 150mm, curved	SIBD9612	426052621TABD01X6	Ila	7
Anatomical tweezers, fine, 145mm	SIBD2166	426052621TABD01X6	Ila	7
Anatomical tweezers, fine tip, 115mm	SIBD2100	426052621TABD01X6	Ila	7
Anatomical tweezers, straight, 180mm	TABD0677	426052621TABD01X6	Ila	7
Anatomical tweezers, medium width, 145mm	TABD0277	426052621TABD01X6	Ila	7
Surgical tweezers, 1x2 teeth, fine,n. SEMKEN, 145mm,	SIBD6699	426052621TABD01X6	Ila	7
Surgical tweezers, strong, 145mm	SIBD6000	426052621TABD01X6	Ila	7
Surgical tweezers, medium width 145mm	TABD5377	426052621TABD01X6	Ila	7
Eyelash tweezers according to Barraquer, octagonal, length 110mm	SIBD4114	426052621TABD01X6	Ila	7
DeBaKey atraumatic tweezers straight 150 mm	TAFB4111	426052621TABD01X6	Ila	7
Jansen tweezers, bayonet-shaped, curved, 160mm	SIBD9613	426052621TABD01X6	Ila	7
Micro Adson tweezers, anatomical, 120mm	TABD2200	426052621TABD01X6	Ila	7
Adson tweezers, anatomical, 120mm	TABD2222	426052621TABD01X6	Ila	7
Micro Adson tweezers, surgical, 120mm	TABD5100	426052621TABD01X6	Ila	7
Adson tweezers, surgical, 120mm	TABD5111	426052621TABD01X6	Ila	7
Anatomical tweezers according to Cushing, 200mm	TABD2399	426052621TABD01X6	Ila	7
Anatomical tweezers, extra wide, 145 mm	SIBD0777	426052621TABD01X6	Ila	7
Surgical tweezers 1x2 teeth 200 mm	TABD6711	426052621TABD01X6	Ila	7
Splinter tweezers straight 90 mm	TABD9644	426052621TABD01X6	Ila	7
Splinter tweezers, straight, 80mm	SIBD9522	426052621TABD01X6	Ila	7

Standard tweezers anatomical straight 130 mm	TABD0455	426052621TABD01X6	Ila	7
Standard tweezers anatomical straight 145 mm	TABD0477	426052621TABD01X6	Ila	7
Anatomical tweezers, 180mm	TABD0488	426052621TABD01X6	Ila	7
Standard surgical tweezers 1x2 teeth 145 mm	TABD5577	426052621TABD01X6	Ila	7

Warnings/Precautions

If you are sensitive to ethylene oxide or stainless steel, consult the medical staff before using the product.

The product must not be used on the central circulatory system and the central nervous system.

The tweezers can only be used for their intended purpose in medical disciplines by trained, medically qualified staff.

Improper storage or handling that deviates from the intended purpose may result in injury to the patient.

Store at **room temperature** and in a place protected from moisture, heat and direct sunlight, otherwise the sterility of the product cannot be guaranteed.



Always keep the tweezers out of the reach of children and other unauthorised persons.

The specified usage times must always be observed. Expired products must be disposed of.

Excessive mechanical forces from the user may cause damage to the product, which may result in injury to the patient.

Moving parts and their unintentional movement can cause injury to the user.

Reusing the tweezers can lead to bio-contamination and an increase in the bacterial load. The same applies if the instruments are used despite the packaging being damaged.



This is a gripping tool; if handled improperly or if the sharp ends are touched, the user may cut themselves or otherwise injure themselves.

Improper disposal of instruments leads to environmental pollution. Please refer to the section on proper disposal.

Intended use

Tweezers are important tools with great benefits that have many applications in the health and medical fields. The tweezers are used to hold, clamp and grasp tissues and materials. They are also

used to hold and pull threads, bandages, swabs; and to hold/remove foreign bodies. They grasp the tissue or material that is to be cut, prepared or sewn, and are accordingly profiled on the inside of the mouth. Tweezers consist of two halves that are pressed together when holding the tissue or material and spring back to their original position immediately when the pressure of the fingers is released.

Due to the numerous possible applications, many different forms of tweezers are available. The respective versions of the products can be seen in the product description in the product list.

Indication: Use in routine surgical procedures in hospitals or medical practices.

Contraindication: The use is contraindicated in cases of general inoperability, previously mentioned intolerances or applications to the central circulatory and nervous systems.

Target patient group: All patient groups can be treated with sterile tweezers

Intended users: The tweezers are intended for use only by trained, medically qualified staff.

Features

There is a sterile barrier for the tweezers.

Correct application

The tweezers can only be used for their intended purpose in medical disciplines by trained, medically qualified staff. The tweezers are only used for a short time, i.e. for a period of up to several hours.

Sterility

If the sterile packaging is damaged, the tweezers must no longer be used. The sterile packaging can only be opened immediately before use. If the sterile packaging is accidentally opened or damaged, please use another product with intact sterile packaging.

The product is sterilised using the ethylene oxide process.



Reuse of the product

The tweezers are designed for single use. Re-sterilisation is prohibited because there is no validated procedure for this process.

Disposal

When disposing of potentially contaminated products, regulations for preventing infection and dealing with microbiological hazards must be observed. This particularly applies to tweezers that have come into contact with tissue or body fluids.

For disposal, country-specific legislation and regulations must be observed. Our instruments are licensed for disposal via the Scholz Recycling System (SReS©).

Version: 3.0

Notice: All serious incidents related to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is located