

Instructions for use for medical tube clamps from Scholz Medical GmbH

Manufacturer

Scholz Medical GmbH
Robert-Bosch-Strasse 1
83052 Bruckmühl
Germany
<https://www.scholz-instruments.com/>
info@scholz-medical.de
+49 8062 8070801



This instruction manual applies to the following products

Description	Product no.	Basic UDI-DI	Category	Rule
Tube clamp, smooth mouth, 180mm, without prongs	TAAN9111	426052621TAAN01YH	Is	5
Tube clamp grooved jaw 180mm	TAAN9122	426052621TAAN01YH	Is	5
Tube clamp grooved jaw 200mm	TAAN9133	426052621TAAN01YH	Is	5
Tube clamp smooth jaw 165mm	TAAN9099	426052621TAAN01YH	Is	5
Tube clamp, strong model, 145mm	TAAN9088	426052621TAAN01YH	Is	5
Tube clamp, plastic, 120mm	TAAN9011	426052621TAAN01YH	Is	5

Warnings/Precautions

If you are sensitive to ethylene oxide or stainless steel, consult the medical staff before using the product.

The product must not be used on the central circulatory system and the central nervous system.

The tube clamps can only be used for their intended purpose in medical disciplines by trained, medically qualified staff.

Improper storage or handling that deviates from the intended purpose may result in injury to the patient.

Store at **room temperature** and in a place protected from moisture, heat and direct sunlight, otherwise the sterility of the product cannot be guaranteed.



Always keep the tube clamps out of the reach of children and other unauthorised persons.

The specified usage times must always be observed. Expired products must be disposed of.

Excessive mechanical forces from the user may cause damage to the product, which may result in injury to the patient.

Moving parts and their unintentional movement can cause injury to the user.

Reusing the tube clamps can lead to bio-contamination and an increase in the bacterial load. The same applies if the instruments are used despite the packaging being damaged.



This is a clamping tool; if handled improperly or the clamping jaws are touched, the user may be crushed or otherwise injured.

Improper disposal of instruments leads to environmental pollution. Please refer to the section on proper disposal.

Intended use

Tube clamps are important tools with great benefits that have many applications in the health and medical fields. The tube clamps are used to clamp various tubes made of medical silicone in the medical and laboratory sectors.

Depending on the tube material, the jaw part must be selected; this is either smooth or cross-grooved. This allows the tube to be gripped and clamped securely, but without being damaged.

Due to the numerous possible applications, many different forms of tube clamps are available. The respective shape and size of the instrument determine its intended use.

Indication: Use in routine surgical procedures in hospitals or medical practices.

Contraindication: The use is contraindicated in cases of general inoperability, previously mentioned intolerances or applications to the central circulatory and nervous systems.

Target patient group: All patient groups can be treated with sterile tube clamps

Intended users: The tube clamps are intended for use only by trained, medically qualified staff.

Features

There is a sterile barrier for the tube clamps.

Correct application

The tube clamps can only be used for their intended purpose in medical disciplines by trained, medically qualified staff. The tube clamps are only used for a short time, i.e. for a period of up to several hours.

Sterility

If the sterile packaging is damaged, the tube clamp must no longer be used. The sterile packaging can only be opened immediately before use. If the sterile packaging is accidentally opened or damaged, please use another product with intact sterile packaging.

The product is sterilised using the ethylene oxide process.



Reuse of the product

The tube clamps are designed for single use. Re-sterilisation is prohibited because there is no validated procedure for this process.

Disposal

When disposing of potentially contaminated products, regulations for preventing infection and dealing with microbiological hazards must be observed. This particularly applies to tube clamps that have come into contact with tissue or body fluids.
For disposal, country-specific legislation and regulations must be observed. Our instruments are licensed for disposal via the Scholz Recycling System (SReS©).

Date of issue of the instructions for use 20/05/2025

Version: 3.0

Notice: All serious incidents related to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is located