

Instructions for use for medical probes from Scholz Medical GmbH

Manufacturer

Scholz Medical GmbH
Robert-Bosch-Strasse 1
83052 Bruckmühl
Germany
<https://www.scholz-instruments.com/>
info@scholz-medical.de
+49 8062 8070801



This instruction manual applies to the following products

Description	Product no.	Basic UDI-DI	Category	Rule
Button probe 1.0 mm x 100 mm	SIBP1010	426052621SIBP013G	Ila	7
Button probe 1.0 mm x 150 mm	SIBP1012	426052621SIBP013G	Ila	7
Button probe 1.5 mm x 150 mm	SIBP1013	426052621SIBP013G	Ila	7

Warnings/Precautions

If you are sensitive to ethylene oxide or stainless steel, consult the medical staff before using the product.

The product must not be used on the central circulatory system and the central nervous system.

The probes can only be used for their intended purpose in medical disciplines by trained, medically qualified staff.

Improper storage or handling that deviates from the intended purpose may result in injury to the patient.

Store at **room temperature** and in a place protected from moisture, heat and direct sunlight, otherwise the sterility of the product cannot be guaranteed.



Always keep the probes out of the reach of children and other unauthorised persons.

The specified usage times must always be observed. Expired products must be disposed of.

Excessive mechanical forces from the user may cause damage to the product, which may result in injury to the patient.

Moving parts and their unintentional movement can cause injury to the user. The probes must not be bent manually or mechanically. This can lead to metal chipping and injury to the patient.

Reusing the probes can lead to bio-contamination and an increase in the bacterial load. The same applies if the instruments are used despite the packaging being damaged.



This is a probing tool; if handled improperly or if the sharp ends are touched, the user may cut themselves or otherwise injure themselves.

Improper disposal of instruments leads to environmental pollution. Please refer to the section on proper disposal.

Intended use

Probes are important tools with great benefits that have many applications in the health and medical fields.

They are used to feel, palpate, measure and/or expand unevenness, abnormalities and injuries in tissue ducts in order to obtain an idea of the depth, size and extent of the diseased part of the body.

Indication: Use in routine surgical procedures in hospitals or medical practices.

Contraindication: The use is contraindicated in cases of general inoperability, previously mentioned intolerances or applications to the central circulatory and nervous systems.

Target patient group: All patient groups can be treated with sterile probes

Intended users: The probes are intended for use only by trained, medically qualified staff.

Features

There is a sterile barrier for the probes.

Correct application

The probes can only be used for their intended purpose in medical disciplines by trained, medically qualified staff. The probes are only used for a short time, i.e. for a period of up to several hours.

Sterility

If the sterile packaging is damaged, the probe must no longer be used. The sterile packaging can only be opened immediately before use. If the sterile packaging is accidentally opened or damaged, please use another product with intact sterile packaging.

The product is sterilised using the ethylene oxide process.



Reuse of the product

The probes are designed for single use. Re-sterilisation is prohibited because there is no validated procedure for this process.

Disposal

When disposing of potentially contaminated products, regulations for preventing infection and dealing with microbiological hazards must be observed. This particularly applies to probes that have come into contact with tissue or body fluids.
For disposal, country-specific legislation and regulations must be observed. Our instruments are licensed for disposal via the Scholz Recycling System (SReS©).

Date of issue of the instructions for use 20/05/2025

Version: 3.0

Notice: All serious incidents related to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is located