clevederm

**EN - Instruction for use**

Product: **clevederm** – Transparent Film Dressing High MVTR

**Remark: The product shape and size can be customized according to applicable requirements.**

**Composition:**

Transparent Film Dressing (High MVTR) is composed of carrier film or paper, PU film, adhesive, release paper and Tape.

**Intended use:**

Transparent Film Dressing is used to fixate the area punctured by injection needles and provide protection.

**Directions for use:**

- Clean and thoroughly dry the application area in accordance with normal procedures.
- Remove the securing strips from the carrier paper and apply the product as required. Ensure that the strips can be used on the catheter wings of hub for improved security. To secure tubing or to stabilize catheter lumens. Remove the carrier from the dressing.
- Remove the release paper and centre the dressing over the catheter insertion site.
- Place the dressing on the catheter hub, adhere the dressing to the patient’s skin and ensure the minimal creasing of the product. Remove the remaining protector paper.
- Remove the transparent carrier by lifting the non-adhesive tab at the end of the dressing and pulling diagonally until the carrier is removed.
- The film around the catheter is molded to complete the application of the dressing.
- Information may be recorded with a ballpoint pen on the documentation label before removing from the carrier and being placed in the patient’s notes, on the skin or on the dressing (without obscuring the view of the catheter insertion site).
- If using more than one dressing, the dressings are not allowed to overlap.
- Remove the dressing by loosening the film and stretching the dressing parallel to the patient’s skin while stabilizing the patient’s skin and catheter with the other hand. Gently and slowly remove the dressing at an angle of 180°.

**Precautions:**

- Pay attention to the applicable scope and use the product appropriately.
- This product is for single use and sterilized by using ethylene oxide.
- For single use only, re-use is prohibited. Re-use may cause cross contamination.
- Check the package before use. Do not use the product if there is any damage or moist.
- Check the expiry date before use.
- Re-sterilization is prohibited.
- Product after use shall be disposed in compliance with local regulations.
- This product should be kept dry and away from sunlight.

**Contraindications:**

Unknown for contraindications.

**Sterilization:**

Ethylene oxide

**Shelf life:**

The product shelf life is 5 years. See the labels of the package for the product lot number and shelf life.

**Storage:**

The product should be kept dry and away from sunlight.

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Shanghai ISO Medical Products Co., Ltd.
Lotus NL B.V.