EN - Instruction for use
Product: **clevefoam** – foam dressing
clevefoam border – foam dressing border

We recommend reading the instructions for use carefully always before using a package of clevocare, as the information is updated regularly.
Intended use:

clevefoam - foam dressing is intended for different phases of wound, especially to highly exudative wounds, (e.g., burns, scalds), cosmetic operation, skin-donor sites.

Indications for Use:
clevefoam - foam dressings are indicated for acute and chronic wounds with moderate to high exudate, including:

Pressure ulcers; Diabetic foot ulcers; Leg ulcers; Post-operative wounds; Skin abrasions; Burns (first and second degree); Donor sites; Traumatic wounds; Shallow, granulating wounds;
Counter-Indications:
Not identified.

Side-effects:

Not identified.

Warning and precautions:
1. Foam dressing is disposable. Do not use if the inner package is damaged.
2. The products have a shelf-life of 3 years. Do not use beyond the expiration date.

3. Foam Dressing is sterilized. Use immediately after opening the packaging, retain a single piece of the product for no longer than 7 days, and no more than 30 days for cumulative use.
4. Foam dressing should be stored in an environment with suitable temperature, protection from light, and ventilation.
5. Immediate medical consultation is required under the following conditions:
- Uncertainty in diagnosis – inability to identify the wound type or underlying medical conditions;
 - Laypersons must not attempt self-application without clinical guidance.
6. Foam dressings are not suitable for arterial bleeds
7. Do not overlay or overlap foam dressings.

Application instruction:

1. **clevefoam - foam dressing non-adhesive:**
1. Cleanse the wound in accordance with local clinical procedures.
 2. Select a suitable size of foam dressing so that the dressing overlaps the wound margins by approximately 2cm. Cut the dressing to size if necessary.
 3. Place the foam side directly to the wound surface.
 4. Secure the dressing with an appropriate bandage or with tape. The tape should only be applied to the edge of the dressing.
2. **clevefoam border - foam dressing adhesive:**
1. Cleanse the wound according to clinical practice. Dry the surrounding skin thoroughly.
 2. Select an appropriate size/shape. The wound pad should cover the dry surrounding skin by approximately 2cm.
 3. Remove the first release film and apply the adherent side to the wound.
 4. Remove the remaining release film(s) and smooth down the border on the skin. Do not stretch the dressing.

The dressing change interval may be several days. Change the dressing before it is fully saturated, at signs of leakage, or as indicated by clinical practice.

Frequency of Dressing Change:

1. The use of foam dressings should always be re-evaluated if wound exudate levels change.
2. During the early stages of wound management, observe the wound daily for signs of infection or other complications (or as recommended by local clinical protocols).
3. Infected wounds should be inspected and treated as per local clinical protocol.
4. The frequency of dressing change is determined by the condition of the wound and local clinical protocols.
5. Dressings can be left in place, undisturbed, for up to 7 days, or until exudate is visible and approaches 2cm from the edge of the dressing, whichever is sooner.

Dressing Removal:
Careful and correct dressing removal is vital to minimize the risk of damaging the wound bed and new granulating tissue;

1. Remove the retention dressing or bandage and carefully lift the dressing away from the wound.
2. The dressing should be removed slowly and carefully. The adhesive foam dressing should be removed following the direction of hair growth.
3. As the dressing is removed, the newly exposed skin should be supported.
4. Soiled dressings should be discarded in the appropriate medical waste receptacle.

Shelf-life and storage information:
Shelf-life: The products have a shelf-life of 3 years.
Storage information: Keep dry. Keep away from sunlight and/or heat. The foam may change color to more yellow when exposed to light, air and/or heat. This has no influence on product properties.
Disposal:
Disposal should be handled according to local environmental procedures.

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