EN - Instruction for use
Product:
clevesil – Silicone interface poliurethane foam dressing
clevesil border – Silicone interface poliurethane foam dressing with border
We recommend reading the instructions for use carefully always before using a package of clevecare, as the information is updated regularly.
Intended use:

Silicone foam dressing is suitable for many types of exuding wounds including leg and pressure ulcers superficial and partial thickness burns, donor sites, postoperative wounds and skin abrasions. The dressing absorbs exudates and maintains a moist wound-healing environment whilst minimizing the risk of maceration.

Indications for Use:
Chronic and acute wounds with light to moderate levels of exudate: Skin abrasions; Pre-tibial lacerations; Pressure ulcers; Diabetic foot ulcers; Leg ulcers; Post-operative wounds; Superficial burns; Donor sites; Traumatic wounds; Shallow, granulating wounds;
Counter-Indications:
Not identified.

Side-effects:
Not identified.

Warning and precautions:

- Silicone Foam dressing is disposable. Do not use if the inner package is damaged.
- The products have a shelf-life of 3 years. Do not use beyond the expiration date.
- This product is sterilized. Use immediately after opening the packaging, retain a single piece of the product for no longer than 7 days, and no more than 30 days for cumulative use.
- Silicone Foam dressing should be stored in an environment with suitable temperature, protection from light, and ventilation.
- Immediate medical consultation is required under the following conditions:
 - Uncertainty in diagnosis – Inability to identify the wound type or underlying medical conditions;
 - Laypersons must not attempt self-application without clinical guidance.
- Silicone Foam dressings are not suitable for arterial bleeds.
- Do not overlay or overlap foam dressings.
- Do not stretch the skin or dressings during application.
- Silicone is hydrophobic and will not adhere to wet skin.

- Application instruction:**
- Cleanse the wound according to clinical practice.
 - Select an appropriate size/shape. The wound pad should cover the dry surrounding skin by approximately 2cm.
 - Remove the first release film and apply the adherent side to the wound.
 - Remove the remaining release film(s) and smooth down the border on the skin. Do not stretch the dressing.
 - Smooth the dressing border down to ensure good adherence with the peri-wound skin.

- Frequency of Dressing Change:**
- The use of silicone foam dressings should always be re-evaluated if wound exudate levels change.
 - During the early stages of wound management, observe the wound daily for signs of infection or other complications (or as recommended by local clinical protocols).
 - Infected wounds should be inspected and treated as per local clinical protocol.
 - The frequency of dressing change is determined by the condition of the wound and local clinical protocols.
 - Dressings can be left in place for up to 7 days, or until exudate becomes visible and approaches 2cm of the dressing margins, whichever occurs first.

- Dressing Removal:**
Careful and correct dressing removal is vital to minimize the risk of damaging the wound bed and new granulating tissue;
- Remove the retention dressing or bandage and carefully lift the dressing away from the wound.
 - The dressing should be removed slowly and carefully following the direction of hair growth.
 - As the dressing is removed, the newly exposed skin should be supported.
 - Soiled dressings should be discarded in the appropriate clinical waste receptacle.

Shelf-life and storage information:
Shelf-life: The products have a shelf-life of 3 years.
Storage information: The product should be stored in a horizontal position. Keep dry. Keep away from sunlight and/or heat. The foam may change colour to more yellow when exposed to light, air and/or heat. This has no influence on product properties.
Disposal:
Disposal should be handled according to local environmental procedures.

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Zhende Medical Co., Ltd.
Shanghai International Holding Corp. GmbH (Europe)